


# ARROWPACK




## ARROWPACK <sup>®</sup>

Tyvek<sup>®</sup> sterilization reels.

Arrowpack Plasma see trough fibre free sterilization packaging reels are only made of the best adapted medical grade materials. For the manufacturing of Arrowpack reels for plasma sterilization, we only use medical grade Tyvek in combination with high quality peelable ecologically produced film.

All Arrowpack reels are produced under cleanroom conditions.



Arrowpack Plasma is designed for hydrogen peroxide sterilization procedure, but can also be used in ETO or other low temperature sterilization procedures.

When printing the indicators, Arrowpack only uses a special non toxic heavy metal free ink .

Arrowpack indicators are printed in between sealing lines in order to avoid any decrease of sealing strength and to prevent migration or contact with the sterile medical devices.

During the sterilization process, Arrowpack indicators change colour from red to blue.

In order to reduce the risk of storage and transport, Arrowpack reels are packed in strong plastic PE covers and high quality protective boxes. Instructions for use are printed symbolically according to standard EN 980.

All Arrowpack sterilization reels are produced according to standard ISO EN11607-1&2 and EN868-5



Dimension	Cat. Nr.	Qty/box
ARROWPACK PLASMA Reel 60mm×70m	55060	12
ARROWPACK PLASMA Reel 75mm×70m	55075	10
ARROWPACK PLASMA Reel 100mm×70m	55100	8
ARROWPACK PLASMA Reel 150mm×70m	55150	6
ARROWPACK PLASMA Reel 200mm×70m	55200	4
ARROWPACK PLASMA Reel 250mm×70m	55250	4
ARROWPACK PLASMA Reel 300mm×70m	55300	2
ARROWPACK PLASMA Reel 350mm×70m	55350	2
ARROWPACK PLASMA Reel 400mm×70m	55400	2
ARROWPACK PLASMA Reel 500mm×70m	55500	2

